

foodproof® *Campylobacter* Detection Kit
- Hybridization Probes (LC 1.x, 2.0, 480 II) -

Order No. R 310 05

Quick Reference Procedure (for LC 480 II)
Version 2, March 2012

A. LightCycler® 480 System Protocol

The following procedure is optimized for use with the LightCycler® 480 Instrument II. Program the LightCycler® 480 Instrument II before preparing the reaction mixes. Use the following LightCycler® 480 Instrument II PCR-program for the foodproof® *Campylobacter* Detection Kit (for details on how to program the experimental protocol, see the LightCycler® 480 System Operator's Manual):

Set-Up					
Detection Format	Block Type	Reaction Volume			
Multi Color HybProbe	96	20 µl			
Filter Combination	dynamic mode, Fluos (465-510), Red 640 (498-640) and Cy 5 / Cy 5.5 (498-660)				
Programs					
Program Name	Cycles	Analysis Mode			
Pre-Incubation	1	None			
Amplification	45	Quantification			
Melting	1	Melting Curves			
Cooling	1	None			
Temperature Targets					
	Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)
Pre-incubation					
Segment 1	37	None	00:02:00	4.4	
Segment 2	95	None	00:10:00	4.4	
Amplification					
Segment 1	95	None	00:00:05	4.4	
Segment 2	59	Single	00:00:35	2.2	
Segment 3	72	None	00:00:15	4.4	
Melting					
Segment 1	95	None	00:00:01	4.4	
Segment 2	40	None	00:00:45	2.2	
Segment 3	75	Continuous		0.19	1
Cooling					
Segment 1	40	None	00:00:30	2.2	

Continued next page



B. Preparation of the PCR Mix

Proceed as described below to prepare a 20 µl standard reaction.

Note: The kit contains two different internal amplification controls which are used depending on the instrument utilized. For the LightCycler® 480 Instrument II use the **foodproof** *Campylobacter* Internal Control (LC 480) (vial 6, black cap). Always wear gloves when handling the multi-well plates. Do not touch the upper surface of the PCR multi-well plate.

1. Use a suitable multi-well plate for LightCycler 480 Instrument II.
2. Thaw the solutions and, for maximal recovery of contents, briefly spin vials in a microcentrifuge before opening. Mix carefully but thoroughly by pipetting up and down.
3. In a 1.5 ml reaction tube, prepare the PCR Mix by adding the following components in the order mentioned below, then mix gently but thoroughly by pipetting up and down:

Component	Volume
foodproof <i>Campylobacter</i> Master Mix, (vial 1, yellow cap)	13 µl
foodproof <i>Campylobacter</i> Enzyme Solution, (vial 2, red cap)	1 µl
foodproof <i>Campylobacter</i> Internal Control (LC 480), (vial 6, black cap)	1 µl
Total volume	15 µl

4. Mix carefully but thoroughly by pipetting up and down. Do not vortex.
 - Pipet 15 µl PCR mix into each well of the PCR plate.
 - For the samples of interest, add 5 µl sample DNA to a well.
 - For the negative control, add 5 µl H₂O, PCR-grade (vial 5, colorless cap) to a well.
 - For the positive control, add 5 µl **foodproof** *Campylobacter* Control Template (vial 4, purple cap) to a well.
 - Seal the plate accurately with an optical sealing foil.
5. LightCycler® 480 Instrument II:
 - Place the plate in a swing bucket centrifuge and centrifuge at 1,500 x g for 30 s.
6. Transfer the capillaries or the PCR plate to the LightCycler® 480 Instrument II.
7. Cycle the samples as described above.

***Note:** For further information please refer to: www.bc-diagnostics.com/?cid=1273065297&lang=1

BIOTECON Diagnostics GmbH
Hermannswerder 17,
14473 Potsdam – Germany
Phone +49 (0) 331 2300-200, Fax +49 (0) 331 2300-299
bcd@bc-diagnostics.com, www.bc-diagnostics.com

foodproof is a registered trademark of BIOTECON Diagnostics. LightCycler is a trademark of Roche Diagnostics.

R_3110_05_21-2 (3)