



Real-Time Cycler Profile Compatibility

Real-Time PCR Cycler	Profile	Order-No.-Ending	Comment
Roche			
LightCycler® 480	LP	-1	
LightCycler® 96	LP	-1	
LightCycler® 2.0	-	-	Special consumables: capillaries
Thermo Fisher Scientific / ABI			
PikoReal® 24	LP	-1	
Applied Biosystems® 7500	RP	-2	
Applied Biosystems® 7500 Fast	LP	-1	
QuantStudio™ 3	LP	-1	Also available with RP block!
QuantStudio™ 5	LP	-1	Also available with RP block!
StepOne™	LP	-1	
StepOnePlus™	LP	-1	
Stratagene / Agilent			
AriaMx	LP / DP	-1 / -3	
Mx3000p	RP	-2	
Mx3005p	RP	-2	
IT-IS			
MyGo Pro	DP	-3	
Dualo 32®	DP	-3	
Dualo 32® R²	DP	-3	
Analytik Jena			
qTower 2.0 / qTower 2.2	RP	-2	
qTower³	LP / RP	-1 / -2	
Bio-Rad			
CFX96™	LP	-1	
CFX96™ Deep Well	LP / RP	-1 / -2	
BMS			
MIC	-	-	Special consumables: 4-strip tubes, clear
Eppendorf			
Realplex 4	RP	-2	
Qiagen (Corbett)			
RotorGene®	-	-	Special consumables: 4-strip tubes, clear and other

Profile: LP = Low Profile, RP = Regular Profile, DP = Different Profile; green highlighted cyclers are available by BIOTECON Diagnostics